

2015	EPOND Price list		
EPOND N°	OEM P/No.	Agilent / Varian Axial Vista, 720-ES and axial-view Liberty	CHF
2050-0710N	2010076500	C-Type ' <b>Caliber</b> ' standard Concentric Glass Nebulizer, 2 ml/min., with 'NebLink' tube (TR-30-C2)	390
2050-0800N		K-Type ' <b>Caliber</b> ' standard Concentric Glass Nebulizer, 2 ml/min., with 'NebLink' tube (TR-50-K2)	390
2050-0710K	2010106800	K-style (conical) <b>Caliber</b> concentric glass nebulizer, (TR-30-K2)	390
2050-0710NT		<b>Tyfoon</b> High-performance concentric glass nebulizer, with NebLink sample tube and argon tube	590
2050-KP50	KP-50-00	<b>Meinhard KR+ Concentric Glass Nebulizer, with PFA argon and sample connectors</b>	390
2050—TQP50	TQP-50-00	<b>Meinhard TQ+</b> Concentric Quartz Nebulizer, with PFA argon and sample connectors	390
2050-M030-C1	TR-30-C1	' <b>Meinhard</b> ' Concentric glass nebulizer, Type C, with PFA argon and sample line connectors	390
2050-M030-A3	TR-30-A3	' <b>Meinhard</b> ' Concentric glass nebulizer, Type A, with PFA argon and sample line connectors	390
2050-M030-K	TR-30-K2	' <b>Meinhard</b> ' Concentric glass nebulizer, Type K, with PFA argon and sample line connectors	390
2050-4250		'NB-Class' ( No-Break ) <b>Boreal</b> Concentric Polymer Nebulizer, for low uptakes	870
2550-2002	20012-0002	HF-resistant PFA Concentric Nebulizer, 100 ul/min uptake version	990
2550-2005	2005-0002	HF-resistant PFA Concentric Nebulizer, 400 ul/min uptake version	990
2550-2028	ME4040	<b>PolyPro-ST</b> HF-res concentric nebulizer, with detachable sample uptake tube	590
2550-2024	ME2040	<b>PFA-ST</b> HF-resistant concentric Nebulizer, with detachable sample uptake tube	990
2030-2028	ES-2030-28	<b>PFA-HS</b> HF-res. High-Solids concentric neb, with detachable sample uptake tube	990
2550-2042	ES-2042	Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.	90
2550-2041	ES-2041	Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.	90
3150-0710	TMM3500	Teflon Burgener <b>MiraMist</b> Parallel-Path Nebulizer	990
3150-0720	PMM4000	PEEK Burgener <b>MiraMist</b> Parallel-Path Nebulizer	990
3150-0730	T2100	<b>Teflon T2100</b> 'Work-horse' Parallel-Path Nebulizer	990
2051-0262	7910043700	Baffled Cyclonic Glass spray chamber, with polymer nebulizer adaptor	390
2051-0260	2010081700	Cyclonic Glass spray chamber, glass nebulizer adaptor with 2 O-rings, ML180021	390
2051-0260-OR		<b>O-rings</b> , packet of 6 for spray chamber neb adaptors	36
2051-0264	ML180022	Baffled Cyclonic Glass spray chamber, glass nebulizer adaptor with 2 O-rings	390
2051-0263	9910070400	<b>Jacketed</b> Cyclonic Glass spray chamber, water-cooled, ML180023	790
2080-0170	9910070400	Varian Water Cooled Glass Spray Chamber	990
2051-026T		HF-resistant Cyclonic spray chamber, PTFE	990
EP-3165-3001		ESI PFA HF-resistant Cyclonic Spray Chamber	1090
EP-2128-0471	ES-2128-0471	ESI 47 mm PFA Spray Chamber and ESI 'PureCap' end cap	1290
EP-2128-1000	ES-2128-1000	Mounting Bracket for PFA Spray Chamber	195
EP-4157-1000	ES-4157-1000	<b>PC3</b> Peltier-cooled Cyclonic Quartz Spray Chamber, with nebulizer and Power Supply	2875
EP-9100-6600		<b>PC3-LT</b> (low-temp.) Peltier-cooled Quartz Spray Chamber, reaches -20 degrees C.	3890
2080-2929		HF-resistant ' <b>Erebus</b> ' Cyclonic spray chamber, with drain tube	990
2070-1000	70-100-0237	Spray Chamber sealing rings, packet of 10	90
<b>2080-0120</b>	<b>2010090400</b>	Varian Vista Axial Quartz Torch with 90 Deg. Bend, Injector ID = 2,3mm, ML180017	350
2080-0125	2010090400-PH	Varian Vista Axial Quartz Torch with 90 Deg. Bend, for organics, Injector ID = 1,4mm (ML180027)	370
2080-0130	2010194800	Varian Vista Axial Low-Flow Quartz Torch with 90 Deg. Bend, High Solids	370
2080-0140	2010081200	Varian, Liberty Axial Quartz Torch 2.3mm Injector	370
		<b>*** discount of 20% for Vista torch orders of 3 or more pcs ***</b>	CHF
2080-0490	9910084700	Varian ICP HF-resistant Axial Torch Kit Includes: Torch Body 2080-0150, Alumina Injector 2080-0450 & Teflon Injector Holder 2080-0480	890

2080-0410	2010073000	Varian Quartz Injector 1.8mm ID	130
2080-0420		Quartz Injector Tube, 0.8mm ID, for organics	195
2080-0150	2010083400	Varian Demountable Liberty Axial Torch, Quartz Body	390
2080-0450	2410067500	Varian Alumina Injector 1.8mm ID, L = 130mm.	290
2080-8801		PTFE Aerosol Transfer Tube, spray chamber to torch, L = 310 mm, OD = 5 mm.	90
EP-1843-0180	ES-1843-0180	ESI Demountable Torch sapphire injector, ID = 1,8mm, O-ring free	990
EP-1842-3020	ES-1842-3020	ESI Demountable Torch platinum injector, ID = 2mm, O-ring free	2290
EP-1842-0000	ES-1842	ESI Demountable Torch Quartz Tubes, for use with O-ring free injector base	490
EP-4399-2300	ES-4399-2300	Male Barbed PFA connector for Transfer Line	130
EP-4399-2002-51	ES-4399-2002-51	PFA Elbow	130
EP-1842-1100		Complete Demountable Torch, with O-ring free injector base and Sapphire Injector	1670
<b>EPOND N°</b>		<b>Desolvators and other accessories</b>	<b>CHF</b>
APEX-Q-VA		<b>APEX-Q</b> Desolvator, with PolyPro concentric nebulizer and <b>XENA</b> uptake meter	POA
EP-4599-3002	ES-4599-300A	<b>Spiro TMD</b> Membrane Desolvation Module ( latest model of ACM ), with Power Supply	3990
8100-2051	PRL-02	<b>'PEARL'</b> Cleaning Device for concentric glass nebulizers – <b>suits all types concentric glass nebs</b>	190
<b>SC-2DX-FAST</b>		<b>SC-2DX Autosampler with 'FAST' sample loop – productivity improvement 2x – 3x</b>	POA
<b>SC-Micro</b>		<b>ESI Compact autosampler with 2 sample racks</b>	POA
<b>oneFAST</b>		<b>'FAST' valve system for productivity gains – adapts to Varian and other autosamplers</b>	7890
V-027-0360		Packet of 5 PFA bottles with PFA caps, 25ml volume, OD=28mm.	120
EP-PERGO	PRG-01	<b>'PERGO' argon humidifier</b>	1099
EPME-PG41-100		100-ml bottle multi-element standard solution for ICP analyses, 28 elements @ 100ppm	360
EPME-PG41-500		500-ml bottle multi-element standard solution for ICP analyses, 28 elements @ 100ppm	790
		<i>reference solutions can be made to measure</i>	
<b>EPOND N°</b>	<b>OEM P/No.</b>	<b>Description, tubes</b>	<b>CHF</b>
<i>See web site</i>		All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12	40
EP-4300-0050		<b>Capillink</b> connectors for PTFE sample capillaries with ID=1,6mm, packet of 6	180
SAM-06A		<b>SamTubes</b> torch side arm connectors, 6mm, pkt 4.	68
SAM-08A		<b>SamTubes</b> torch side arm connectors, 8mm, pkt 3.	78
		<b>Autosamplers and ICP / AAS Standard and Reference Solutions</b>	
SC-8000-XXXX		<i>Complete range of ESI SC-series autosamplers - please ask for details</i>	POA
PVCSD-40	3710034400	Tygon peristaltic-pump tubing, packet of 12, ID = 1,02mm	40
PVCSD-65	3710034600	Tygon peristaltic-pump tubing, packet of 12, ID = 1,65mm	40
Tube-SolF-30A		SolventFlex, 2-Stops 72 mm, 190 mm, 0.030" (0.76 mm)	60
SanSD-60-3		Santoprene 3-Stop Tubing 0.060" Yellow/Blue	95
ICP-standards		<i>Certified ICP / AAS Standard solutions - <a href="http://www.e-pond.biz/standards">http://www.e-pond.biz/standards</a></i>	
		<b>*** custom standards made to your specifications ***</b>	
		<b>www.e-pond.biz</b>	