

2015	Components for Agilent HP4500, HP7500 and HP7700/8800 ICP-MS		
EPOND No.	OEM P/No.	Description	CHF
2050-0710NT		'Tyfoon' high-TDS concentric glass nebulizer, with NebLink sample tube	590
2050-0770N	G1820-65030	EPOND A-type Concentric Glass Nebulizer, with NebLink sample tube	390
2050-0770Q	G1820-65138	EPOND A-type Concentric Quartz Nebulizer, with NebLink sample tube	590
2050—TQP50	TQP-50-00	Meinhard TQ+ Concentric Quartz Nebulizer, with PFA argon and sample connectors	530
2050-M030-C1	TR-30-C1	'Meinhard' Concentric glass nebulizer, Type C3, with PFA argon and sample line connectors	390
2050-M030-A3	TR-30-A3	'Meinhard' Concentric glass nebulizer, Type A3, with PFA argon and sample line connectors	390
2550-2020	ES-2020-0002	HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 20 ul/min uptake version	1090
2550-2000	ES-2000-0002	HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 50 ul/min uptake version	990
2550-2002	ES-2002-0002	HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 100 ul/min uptake version	990
2550-2003	ES-2003-0002	HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 200 ul/min uptake version	990
2550-2005	ES-2005-0002	HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 400 ul/min uptake version	990
2550-2028	ES-4040-00xx	PolyPro-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube	590
2550-2024	ES-2040-00yy	PFA-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube	990
2550-2024-07	ES-2040-7000	PFA-ST-APEX concentric nebulizer for use with APEX desolvator	1050
2030-2028	ES-2030-28	PFA-HS HF-res. High-Solids concentric neb, with detachable sample uptake tube	990
2550-2045	ES-2045	Detachable sample capillary tube, ID = 0,15mm., nominal uptake 20 ul/min.	90
2550-2043	ES-2043	Detachable sample capillary tube, ID = 0,20mm., nominal uptake 50 ul/min.	90
2550-2042	ES-2042	Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.	90
2550-2041	ES-2041	Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.	90
		Capillaries with integrated probes - available for all common autosamplers	
3150-0710	G3161-80001	Burgener PTFE MiraMist Parallel-Path Nebulizer for Agilent ICP-MS, TMM3500	990
3150-0720	G3161-80000	Burgener PEEK MiraMist Parallel-Path Nebulizer for Agilent ICP-MS, PMM4000	990
2052-0300	G1820-65121	PTFE Endcap for CG Nebulizer, 35/6 mm adaptor	295
EP-3262-1111	ES-3262-1111	'SSI' Quartz High-stability Dual spray chamber	1970
EP-2122-5351	ES-2122-5351	PFA 35-mm HF-resistant spray chamber with O-ring-free End-Cap	1090
2090-0100		HP Pyrex Spray Chamber 180 deg.	630
2090-0110	G1820-65336	HP Pyrex Spray Chamber 90 deg.	630
2090-0120		HP Quartz Scott Spray Chamber 180 deg.	690
2090-0130	G1820-65337	HP Quartz Scott Spray Chamber for HP4500 and HP7500	690
2090-0140	G1820-80417	HP Quartz Spray Chamber 90 deg. Ext. Inner Tube	690
2090-0150		HP Pyrex Scott Spray Chamber	630
2090-0130	G3280-80008	Quartz Spray Chamber for Agilent 7700 / 8800, Meinhard ML190060	690
2090-016Q		Quartz Cyclonic Spray Chamber	690
ES-4157-1000		PC3 Peltier-cooled Quartz Spray Chamber, with Power Supply	3290
2090-0050	G1820-65028	HP4500 Standard Quartz Torch with 2.5mm Injector	490
2090-0060	G1833-65423	HP7500 Standard Quartz Torch with 2.5mm Injector, original version	490
2090-0062	G3270-80027	HP7500 Standard Quartz Torch with 2.5mm Injector, HMI compatible, vsn 1	490
2090-0010	G3270-67002	HP7500 Standard Quartz Torch with 2.5mm Injector, HMI compatible / ES-1302-7511	490
1302-4510	EP1302-4510	HP Quartz Torch with 1.5mm Injector ID, for HP4500 (for organics)	390
2090-0070	G1833-65424	HP Quartz Torch with 1.5mm Injector ID, for HP7500 (for organics) / ES-1302-7510	390
2090-0510	G3280-80001	Quartz Torch for HP7700, injector ID = 2,5mm - Meinhard ML190050	390
2090-0540	G3280-80004	Quartz Torch for HP7700, tapered injector ID=1,5mm	390
1302-7700	ES-1302-7700	Quartz Outer Tube for O-ring-free HP7700 torch	340
1303-7180	ES-1303-7180	1,8mm Sapphire Injector, O-ring free, for HP7700/8800	650
ML 190065	G3270-80024	HP770 Quartz T-connector with optional gas inlet, ball to socket, L=103mm	85
2090-0080		HP Demountable Torch, Quartz Body	340
2090-0200		HP Alumina Injector, for 2090-0080	270
2090-020Q		HP Quartz Injector, for 2090-0080	190
2090-0250	G1833-65421	Quartz Bonnet for HP7500 and HP7700	140
2090-0300		HP Teflon Injector Holder	195
1302-4545		O-ring-free Demountable Torch for HP4500, with quartz injector in PFA base	790
1302-4500	ES-1302-4500	Quartz Outer Tubes for O-ring-free HP4500 torch	390
1302-7575		O-ring-free Demountable Torch for HP7500, with quartz injector in PFA base	790
1302-7500	ES-1302-7500	Quartz Outer Tubes for O-ring-free HP7500 torch	390
1323-0000	ES-1323-0000	Quartz Injector for O-ring-free Agilent torch, ID = 1,0 / 1,5 / 2,0 or 2,5 mm.	390
1303-0000	ES-1303-0000	Sapphire Injector for O-ring-free Agilent torch, ID = 1,0 / 1,5 / 1,6 or 1,8 mm.	790
1313-0000	ES-1313-0000	Platinum Injector for O-ring-free Agilent torch, ID = 1,5 / 1,8 / 2,0 or 2,5 mm.	2190
2090-0320		HP Quartz Socket Adapter	190

2090-0360		HP Quartz Ball Right-angle Adaptor	190
2090-0380	G1820-80237	HP Quartz Ball & Socket Straight Adaptor, for spray chamber to torch	190
MC-65419-Ag	G1833-65419	Silver (Ag) Shield Plate for quartz torch, for HP7700 ICP-MS	355
MC-65357	G1820-65357	Platinum (Pt) Shield Plate for quartz torch, for HP4500 and HP7500 ICP-MS	810
MC-65419	G1833-65419	Platinum (Pt) Shield Plate for quartz torch, for HP7700 ICP-MS	810
EPOND N°		Other sample-introduction systems and accessories	CHF
ACM-E2		ESI Nafion Active Membrane Device , with DC Power Supply	4890
EP-4599-3002	ES-4599-300A	Spiro TMD Membrane Desolvation Module (latest model of ACM), with Power Supply	4995
SC-1207DX		Enclosure for SC-2DX autosampler	975
8100-2051	PRL-02	'PEARL' Cleaning Device for concentric glass nebulizers	215
EPMS-PG48-100	8500 – 6943	100 ml. bottle - 10 mg/l. Tuning Solution : Ce, Ti, Y, Li in 2% HNO3	195
EPMS-PG48-250		250 ml. bottle - 10 mg/l. Tuning Solution : Ce, Ti, Y, Li in 2% HNO3	270
EP2768.D001.2N.L5	5185 – 5959	500 ml. Bottle - 1 ug/l. Tuning Solution : Ce, Co, Li, Mg, Ti, Y in 2 % w/w HNO3	198
		reference solutions made to measure - on request	
EPOND No.	Agilent No.	Description of Cones	CHF
HP1001-Ni	G1820-65238	Nickel Sampler / ES-3000-1601	495
HP-1001T-Ni	G1820-65480	Nickel Sampler, T Mode / ES-3000-1602	550
HP1002-Ni	G1820-65265	Nickel Skimmer, for HP 4500 / ES-3000-1603	405
HP1002-Ni	G1820-65050	Nickel Skimmer, for HP 7500 / ES-3000-1610	405
HP1002C-Ni	G1833-65497	Nickel Skimmer, for HP 7500c / ES-3000-1604	405
HP1002CE-Ni	G3270-65024	Nickel Skimmer, for HP 7500ce / ES-3000-1611	395
HP-1002T-Ni	G1820-65481	Nickel Skimmer, T Mode / ES-3000-1605	455
HP-5001	G1820-65025	O-ring for sampler cone / ES-3000-1653	10
HP1006-Pt	G1820-65239	Pt Sampler cone with 10mm insert (4500 & 7500)	2095
HP1006A-Pt	G1820-65360	Pt Sampler cone with 18mm insert (4500 & 7500)	5790
HP1008-Pt	G1820-65237	Pt Skimmer, 7.2mm insert	1455
HP1008C-Pt	G1833-65092	Platinum Skimmer for HP 7500c, collision-cell ICP-MS	1340
HP1008CS-Pt	G1883-65132	Pt Skimmer, for HP 7500cs and 7500ce	1340
HP2001-Ni	G3280-67040	Nickel Sampler cone for 7700	585
HP2000P-Ni	G3280-67061	Nickel-plated Sampler cone for 7700	670
HP2002S-Cu	G3280-67067	Copper Skimmer cone for 7700s	535
HP2002S-Ni	G3280-67066	Nickel Skimmer cone for 7700s	475
HP2002X-Ni	G3280-67041	Nickel Skimmer cone for 7700x	475
HP2006-Pt	G3280-67036	Platinum Sampler cone for 7700	2255
HP2006A-Pt	G3280-67056	Platinum Sampler cone for 7700, with 18mm Pt insert	6150
HP2008S-Pt/Cu	G3280-67064	Platinum Skimmer cone for 7700s, with Cu base	1455
HP2008S-Pt/Ni	G3280-67065	Platinum Skimmer cone for 7700s, with Nickel base	1590
HP2008X-Pt	G3280-67060	Platinum Skimmer cone for 7700x	1440
HP2008X-Pt/Ni	G3280-67063	Platinum Skimmer cone for 7700x, Nickel base	1440
G3280-20504	G3280-20504	Retaining Ring for Agilent 7700 sampler cone	275
HP6001	G3280-67009	Gaskets for Agilent 7700 sampler cone, packet of 3	62
2085-2206-05		Diamond paste for polishing and cleaning cones, 3-micron, 5-ml syringe	51
2085-2206-18		Diamond paste for polishing and cleaning cones, 3-micron, 18-ml syringe	93
EPOND No.	OEM No.	Description of ICP-MS detector	CHF
EP14753	AF-0573	HP4500 Electron Multiplier Detector, Agilent P/No. G1820-80386	2290
EP14222	AF-0222	HP7500 Electron Multiplier Detector, Agilent P/No. G1833-65420	2290
EP14HPK	G1833-65575	HP7500 Electron Multiplier Detector, Agilent P/No. G1833-65575	2090
EPOND N°	OEM P/No.	Description, peristaltic pumps and tubing	CHF
ES-4398-1076	Tube-PVCSD-30	Peristaltic-Pump 2-tag tubes, pkt 12, ID=0.76mm, EPOND P/No. EPT-0008	36
ES-4398-1120	Tube-PVCSD-90	Peristaltic-Pump 2-tag tubes, pkt 12, ID=2.29mm, EPOND P/No. EPT-0022	36
EPT-0059	G1833-65569	Agilent peristaltic-pump tubes, 3 tags, ID = 1,02 mm	48
EP-PERGO		'Pergo' argon humidifier - keeps nebulizers running smoothly : 10% discount when bought with nebulizer	1099
ERT-040		EPOND 'RubeTubes', sample tube connectors, packet of 40	165
ICP / AAS Standard, Reference and Custom-made Solutions			
ICP-standards		Complete range of Certified Mono-Element and Multi-Element standard solutions	POA
www.e-pond.biz			