

EPOND S.A. Price List, PE Optima Radial View ICP			
2015			
EPOND N°	OEM P/No.	Description	CHF
2050-0700N	N810-2011	Standard Concentric Glass Nebulizer, with <i>NebLink</i> sample tube	390
2050-0700NQ		Standard Concentric Quartz Nebulizer, with ' <i>NebLink</i> ' sample tube	670
2050-0710N	0047-2022	Hi-flow Concentric Glass Nebulizer, with <i>NebLink</i>	390
2050-0710N	N068-1574	Standard Concentric Glass Nebulizer, with <i>NebLink</i> sample tube	390
TR-30-C1	N8102011	' <i>Meinhard</i> ' Concentric glass nebulizer, Type C	390
TR-30-A3	WE02-4371	' <i>Meinhard</i> ' Concentric glass nebulizer, Type A	390
TR-30-K2	N068-1574	' <i>Meinhard</i> ' Concentric glass nebulizer, Type K	390
2050-0710NT		High-Performance ' <i>Tyfoon</i> ' Concentric Glass Nebulizer, with <i>NebLink</i>	590
2550-2005		HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 360 ul/min uptake version	990
2550-2028	ME4040	<i>PolyPro-ST</i> HF-resistant concentric Nebulizer, with detachable sample uptake tube	590
2550-2024	ME2040	<i>PFA-ST</i> HF-resistant concentric Nebulizer, with detachable sample uptake tube	990
2550-2042	ES-2042	Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.	90
2550-2041	ES-2041	Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.	90
3150-0710	TMM3500	Burgener Teflon <i>MiraMist</i> Parallel-Path Nebulizer	990
3150-0720	PMM4000	Burgener PEEK <i>MiraMist</i> Parallel-Path Nebulizer	990
3150-0730	T2100	Teflon T2100 'Work-horse' Parallel-Path Nebulizer	970
2052-0300		PTFE Endcap for CG Nebulizer, 35/6 mm adaptor	240
2030-0150		Optima 3000 Radial Spray Chamber, Top Port Angled	490
2030-0160		Scott (Coaxial) Spray Chamber, Optima 3000 Radial	490
2030-0190		Cyclonic Baffled Spray Chamber for Optima 3000 radial	490
2030-0200		MLN Cyclonic Spray Chamber, Optima 3000 Radial	490
2030-0194		' <i>Libra</i> ' baffled cyclonic glass spray chamber, 50 ml. With 11-mm ID output tube	490
2031-0350	6808-0204	Spray-Chamber Torch adaptor for cyclonic spray chambers	330
EP-4157-1000		<i>PC3</i> Peltier-cooled <i>Thermostatted</i> Cyclonic Spray Chamber, with 12-Volt Power Supply	2470
EP-9100-5700		<i>PC3-LT</i> (low-temp.) Peltier-cooled Quartz Spray Chamber, reaches -20 degrees C.	3390
2031-0040	N069-0568	Optima 3000 Radial Quartz Torch	390
2031-004W		Optima 3000 Radial Torch, for Wear Metals application	390
2031-0050	N069-5295	Torch Bonnet for Optima 3000 radial, before sept. 1994	95
2031-0130	N069-5456	Optima 3000 Radial Qtz Torch Bonnet - post Sept. 1994	95
2032-0330	N068-1631	Qtz Sample Injector Optima 3000 Radial, ID=1.2mm	130
2032-0430	N069-5448	Qtz Sample Injector Optima 3000 Radial, ID=1.6mm	130
2032-0400	N069-5441	Qtz Sample Injector Optima 3000 Radial, ID=2.0mm	130
2032-0510	N058-2186	Alumina Injector, ID=0.85mm	290
2032-0520	N058-2185	Alumina Injector, ID=1.5mm	290
2032-0530	N058-2184	Alumina Injector, ID=2.0mm	290
EPOND P/N°		Other Sample Introduction Accessories	CHF
APEX-Q-PE		<i>APEX-Q</i> Desolvator, High-Sensitivity Sample Introduction System – raises sensitivity up to 8x	11950
EP-4599-3002	ES-4599-300A	<i>Spiro TMD</i> Membrane Desolvation Module with Power Supply	3990
8100-2051		' <i>PEARL</i> ' Cleaning Device for concentric glass nebulizers	190
2045-0700		' <i>PERGO</i> ' Argon Humidifier, new technology, uses vapour permeation	990
EPOND N°	OEM P/No.	Description, tubes	CHF
EP-MP2-04		<i>MP2</i> Precision peristaltic pump, 4-channels, can be integrated in ICP panel	2390
See web site		All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12	48
ERT-040		EPOND ' <i>RubeTubes</i> ', sample tube connectors, packet of 40	165
ICP-standards		Certified ICP / AAS Standard solutions - <a href="http://www.e-pond.biz/standards">http://www.e-pond.biz/standards</a>	
<a href="http://www.e-pond.biz">www.e-pond.biz</a>			