

2015		EPOND Switzerland Parts List	
EPOND N°	OEM P/No.	Parts for PE-Sciex Elan 6000, 6100, DRC 9000 and NexION	CHF
		<i>note : NB-class nebulizers are HF-resistant</i>	
2050-0710N		Standard ' Caliber ' Concentric Glass Nebulizer, with 'NebLink' sample tube	390
TR-30-C3	N8102011	' Meinhard ' Concentric glass nebulizer, Type C3	390
TR-30-A3	WE02-4371	' Meinhard ' Concentric glass nebulizer, Type A3	390
TR-30-K2	N068-1574	' Meinhard ' Concentric glass nebulizer, Type K	390
2050-0710NT		Tyfoon high-TDS concentric glass nebulizer, with 'NebLink' sample tube and argon connector	490
TR-50-C0.5	N8145012	Low-uptake concentric glass nebulizer for PE NexION or PE Elan	410
TQ-50-A0.5	N8145011	Low-uptake concentric quartz nebulizer for PE NexION or PE Elan	440
TQP-30-A3	WE024371	Meinhard A-type Quartz Concentric nebulizer, with 'NebLink' sample tube and argon connector	570
2050-0750NQ		Concentric Quartz Nebulizer, with 'NebLink' sample tube and argon connector	570
ME2040-7000	N0777375	Apex-ST PFA MicroFlow Nebulizer	1200
ME2002-7000	N0777376	High temperature APEX-100 PFA Nebulizer	1200
2550-2020	N8122350	HF-resistant MicroFlow PFA Concentric Nebulizer, 20 <i>ul/min</i> uptake version	990
2550-2000	ES-2000-0002	HF-resistant MicroFlow PFA Concentric Nebulizer, 50 <i>ul/min</i> uptake version	990
2550-2002	N8122351	HF-resistant MicroFlow PFA Concentric Nebulizer, 100 <i>ul/min</i> uptake version	990
2550-2003	ES-2003-0002	HF-resistant MicroFlow PFA Concentric Nebulizer, 200 <i>ul/min</i> uptake version	990
2550-2005	ES-2005-0002	HF-resistant MicroFlow PFA Concentric Nebulizer, 400 <i>ul/min</i> uptake version	990
2550-2028	ES-4040-28	PolyPro-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube	590
ES-2040-54	N8122192	PFA-ST HF-resistant Concentric PFA Nebulizer, with detachable sample uptake tube	990
2550-2024	ES-2040-27	PFA-ST HF-resistant Concentric PFA Nebulizer, with detachable sample uptake tube, nominal uptake 100uL/min	990
2030-2028	ES-2030-28	PFA-HS HF-res. High-Solids concentric neb, with detachable sample uptake tube	990
2550-2045	ES-2045	Detachable sample capillary tube, ID = 0,15mm., nominal uptake 20 <i>ul/min</i> .	90
2550-2043	ES-2043	Detachable sample capillary tube, ID = 0,20mm., nominal uptake 50 <i>ul/min</i> .	90
2550-2042	N8122383	Detachable sample capillary tube, ID = 0,25mm., uptake 100 <i>ul/min</i> ., ES-2042	90
2550-2041	N8122384	Detachable sample capillary tube, ID = 0,50mm., uptake 400 <i>ul/min</i> ., ES-2041	90
		Capillaries with integrated probes - available for all common autosamplers	
3150-0710	TMM3500	Burgener PTFE MiraMist Parallel-Path Nebulizer	990
3150-0720	PMM4000	Burgener PEEK MiraMist Parallel-Path Nebulizer	990
3150-5000	AM-5000	Burgener AriMist nebulizer	990
2052-0300		PTFE Endcap for CG Nebulizer, 35/6 mm adaptor, similar to N0680504	240
2030-019MS	N812-0150	' Althea ' Cyclonic glass Spray Chamber for Elan 5000/6000/9000/NexION, # 2030-0192G	450
2030-0192G	N812-0150	' Althea ' Cyclonic glass Spray Chamber for Elan 5000/6000/9000/NexION	450
2030-019MSJ		Jacketed ' Althea ' Cyclonic Spray Chamber for Elan 5000/6000/9000/NexION	650
2030-019MSQ		Quartz ' Althea ' Cyclonic Spray Chamber for Elan 5000/6000/9000/NexION	650
2030-019QMSJ		Jacketed Quartz Cyclonic Chamber, Elan 5000/6000/9000/NexION	650
3154-0193		Quartz Mini-Cyclonic quartz spray chamber, 'Helical' nebulizer port, RS12/5 coupling	990
3154-1119	EP-3154-1119	O-ring-free Quartz ' Libra ' Spray Chamber for Elan DRC, with auxiliary gas input	990
2040-0200	WE02-5221	Quartz ' Libra ' Baffled Spray Chamber for Elan 6100DRC	890
2040-020JQ		Jacketed Quartz Cyclonic Spray Chamber, Elan 6100DRC	990
EP-3154-3011	ES-3154-3011	HF-resistant Cyclonic Spray Chamber for Elan 5000/6000	1140
EP-5010-0350	ES-5010-0350	PFA 35/6mm O-ring free ESI Spray Chamber Adaptor, with Aux gas input port	490
EP-4157-100C	N8122410	PC3 Peltier-cooled ESI Quartz Spray Chamber + Power Supply, Cassette Torch Mount	2995
EP-4157-100T	N8122389	PC3 Peltier-cooled ESI Quartz Spray Chamber + Power Supply, Twist Torch Mount	2995
EP-9100-5400	ES-9100-5400	PC3-LT (low-temp.) Peltier-cooled Quartz Spray Chamber, reaches -20 degrees C.	3890
EP-2100-0472	N8122356	47 mm ESI PFA Coaxial (Scott type) Spray Chamber, without nebulizer adapter	1140
EP-2100-9479	ES-2100-5472	47 mm ESI PFA Coaxial (Scott type) Spray Chamber, with Elan torch & neb. adapters	1570
ES-2105-2103	ES-2105-2103	HF-resistant Drain Adapter kit, with Luer connectors	70
ES-2044-0002	ES-2044-0002	HF-resistant Drain Adapter kit, with threaded connector	70
ES-5010-0470	N8122357	PFA End-Cap for 47mm Scott type spray chamber	450
7808-0404		RS13/5 ELAN Elbow Adaptor	165
2040-0100	N812-2006	Elan Standard Quartz Torch, for Elan 5000/6000/6100/9000, ML140010	225
2040-0308	WE02-3948	Tapered Quartz Injector with ball joint 2mm (for Elan 6100 & DRC), ML140030	330
2040-0370	WE02-3951	Ball Joint Injector Support Adapter for Elan DRC, suits #2040-0308	220
2040-0350	N8120016	Injector Support, Torch Base for Elan 5000/6000	400
2040-0360	N812-0100	PE/Sciex O-Ring Adapter Kit	42
2032-0330	N068-1631	Standard PE Quartz Injector without support, ID = 1,2mm for all Elan models	155
2032-0340	0047-3292	Standard PE Quartz Injector without support, ID = 1,6mm for all Elan models	155
2032-0350	N812-5029	Standard PE Quartz Injector without support, ID = 2,0mm for all Elan models	155
MA2186	N8126039	Alumina Injector, 0.85mm	165

MA2185	N8126040	Alumina Injector, 1.5mm	165
2032-0530	N8126041	Alumina Injector, ID = 2,0 mm., MA2184	165
EP-1123-1101	WE027030	Quartz Injector with PFA base, ball-joint, bayonet mount, ID=1,0mm	355
EP-1123-1151	WE027005	Quartz Injector with PFA base, ball-joint, bayonet mount, ID=1,5mm	355
EP-1123-1201	WE023948	Quartz Injector with PFA base, ball-joint, bayonet mount, ID=2,0mm	355
EP-1123-3150	ES-1123-3150	Quartz Injector with ball-joint adaptor, cassette mount, ID=1,5 mm	400
EP-1103-0180	ES-1103-0180	Sapphire Injector, bayonet mount, ID = 1,8 mm (ID of 1,6 and 1,0 mm also available)	930
EP-1113-0100	ES-1113-0100	Platinum Injector, bayonet mount, ID=1,0 mm with PFA support	1915
EP-1113-0150	ES-1113-0150	Platinum Injector, bayonet mount, ID=1,5 mm with PFA support	1915
EP-1113-0201	ES-1113-0200	Platinum Injector, bayonet mount, ID=2,0 mm with PFA support	1915
ES-1101-0107	ES-1101-0107	Platinum shield disc for Perkin Elmer Elan PFA injector base	210
ES-1101-3107	ES-1101-3107	Platinum Shield Disc for ESI O-ring free injectors	210
EP-1103-2180	ES-1103-2180	Sapphire Injector, cassette mount, ID = 1,8 mm (ID of 1,6 and 1,0 mm also available), PE # N8122411	990
EP-1113-3100	ES-1113-31010	Platinum Injector, ID = 1,0 mm with PFA support, cassette mount	1915
EP-1113-3150	ES-1113-3150	Platinum Injector, ID = 1,5 mm with PFA support, cassette mount	1915
EP-1113-3200	ES-1113-3200	Platinum Injector, ID = 2,0 mm with PFA support, cassette mount, PE # N8122412	1915
ES-1113-2200	N8145111	Platinum Injector, for O-Ring-Free (ORF) torch, cassette mount, ID = 2mm	1880
ES-1103-2180	N8145072	Sapphire Injector, for O-Rin-Free (ORF) torch, cassette mount, ID = 1,8 mm (ID of 1,6 and 1,0 mm also available)	780
ES-1124-3100	N0777410	PFA cassette mount quartz injector for NexION , ID=1 mm	450
ES-1124-3150	N8122413	PFA cassette mount quartz injector for NexION , ID=1,5 mm	450
ES-1124-3200	N0777457	PFA cassette mount quartz injector for NexION , ID=2 mm	450
EP-1102-1000	ES-1102-1000	Precision Quartz Torch for O-ring-free injectors (injectors listed below)	490
EP-1124-1200	ES-1124-1200	O-ring free Quartz Injector with PFA support, bayonet mount, ID=2mm	490
EP-1103-1180	ES-1103-1180	O-ring free Sapphire Injector with PFA support, bayonet mount, ID=1,8 mm	930
EP-1113-1200	ES-1113-1200	O-ring free Platinum Injector with PFA support, bayonet mount, ID=2mm	1915
ES-2170-3472	N8122414	PFA Teflon Sample Introduction Kit, Cassette Torch Mount, Sapphire Injector 1,8mm	4995
ES-2171-3472	N8122415	PFA Teflon Sample Introduction Kit, Cassette Torch Mount, Platinum Injector 2,0mm	5490
ES-2170-1472	N8122362	PFA Teflon Sample Introduction Kit, Twist Type Torch Mount, Sapphire Injector 1,8mm	5175
ES-2171-1472	N8122363	PFA Teflon Sample Introduction Kit, Twist Type Torch Mount, Platinum Injector 2,0mm	5490
EPOND N°	OEM P/N°	Probes for PerkinElmer Model AS-9X Autosamplers	CHF
		<i>HF-resistant sampler probes, polyimide tube with PTFE covering</i>	
EP-5007-9150	ES-5007-9150	PE AS-9X series probe with polyimide core. Uptake range 10 – 80 ul/min.	330
EP-5007-9200	ES-5007-9200	PE AS-9X series probe with polyimide core. Uptake range 20 – 160 ul/min.	330
EP-5007-9250	ES-5007-9250	PE AS-9X series probe with polyimide core. Uptake range 50 – 400 ul/min.	330
EP-5007-9500	ES-5007-9500	PE AS-9X series probe with polyimide core. Uptake range 100 – 800 ul/min.	330
EP-5007-9800	ES-5007-9800	PE AS-9X series probe with polyimide core. Uptake range 200 – 1800 ul/min.	330
EP-5009	ES-5009	PerkinElmer Probe Holder : for self-aspirating sample probes	130
EPOND N°	OEM P/N°	Parts for other autosamplers	CHF
SC-5037-3995-150	SC-5037-3995-150	FAST Sample probe with 1mm ID, grey tag, 150 cm length	230
SC-5037-3755-150	SC-5037-3995-150	FAST Sample probe with 0,8mm ID, blue tag, 150 cm length	230
SC-1207-DX-1200	SC-1207-DX-1200	Enclosure for SC-2 Autosampler with ULPA Filter. Includes both enclosure SC-1207 and SC-0602 ULPA (Ultra-Low Particle Arrester) filter unit (99.999% efficient at 0.12 micron particle size).	2480
EPOND N°	OEM P/N°	DESCRIPTION of ICP-MS DETECTOR	CHF
EPM14210	N8125001	Detector for Elan 6000 and 6100DRC	2490
EPM14217	N8125050	Improved Dual Detector for Elan 6100DRC and Elan 9000, post-April, 2005	2540
MC-45000	MC-45000	PE NexION 300 ICP-MS Detector	1730
EPOND N°		Exo-Chambers and other accessories	CHF
ES-2999-5420		RF Load Coil, for Elan 6000, Gold-plated, includes fittings	375
S032110	WE012411	RF Load Coil, Copper, with fittings, for Elan 6000 : PE # WE012411, MC-12411	335
MC-21816	WE021816	RF Load Coil, for Elan 6100, Elan DRC, and NexION : Copper	335
ACM-E2		ESI Nafion <i>Active Membrane Device</i> , with DC Power Supply	4065
EP-4599-3002	ES-4599-300A	Spiro TMD Membrane Desolvation Module (latest model of ACM), with Power Supply	4270
PRG-01		' PERGO ' argon humidifier, suits all nebulizers, NO pressure vessel – based on membrane permeation	995
8100-2051	PRL-02	' PEARL ' Cleaning Device for ALL types of concentric glass nebulizers	215
MS29A0.D01.1N.L25	N8125031	250-ml bottle Elan 6100 / DRC / 9000 Detection Limit Solution : 10 ppb of Be, Co, In and U	195
MS29A0.D01.1N.L5		500-ml. Bottle of Elan Detection Limit Solution : 10 ppb of Be, Co, In and U in 1% HNO3	285
MSB3B8.D01.1N.1L	N8125030	1-Litre bottle of Setup / Stability / Masscal standard solution , 9 elements at 10 ppb concentration each => Ba, Cd, Ce, Cu, In, Pb, Mg, Rh and U, in 1% HNO3	355
8AC1.D01.1N.L5	N8125040	500-ml bottle of SmartTune ICP-MS standard solution for Elan 6000, DRCE, 9 elements in 1% HNO3 => Ba, Be, Ce, Co, In, Mg, Pb, Rh and U	410

0E83.D01.05N.L5	N8125041	500-ml bottle of SmartTune ICP-MS standard solution for Elan DRC, DRC+ and DRC II, 10 elements in 0,5% HNO3 => Ba, Be, Ce, Co, In, Fe, Pb, Mg, Th and U	445
EPOND N°	OEM P/N°	DESCRIPTION of CONES	CHF
SC2011-Ni	WE021140	Elan 6000/6100 and 9000 Nickel Sampler cone, ES-3000-1709, MC-21140	330
ES-3000-1709-S	WE021140	Elan 6000/6100 and 9000 Sampler cone, 100% Nickel (no copper)	330
SC2012-Ni	WE021137	Elan 6000/6100 and 9000 Nickel Skimmer cone, ES-3000-1710, MC-21137	395
MC-33995	W1033995	Al Hyper Skimmer Cone for PE NexION, MC-33995	245
MC-33612	W1033612	Ni Sampler Cone for PE NexION, MC-33612	355
MC-40148	WE1040148	Gasket for PE NexION Sampler Cone, pack of 3	26
MC-26356	W1026356	Nickel Skimmer cone for PE NexION	375
MC-19737	9902123	Viton O-ring for Hyper Skimmer, packet of 2	24
MC-33614	W1033614	Pt sampler Cone, Ni base for Perkin NexION 300s	2790
MC-26907	W1026907	Platinum skimmer cone for NexION	1090
SC2013-Pt	WE027802	Elan 6000/6100 and 9000 Pt Sampler cone, MC-27802	2390
SC2014-Pt	WE027803	Elan 6000/6100 and 9000 Pt Skimmer cone, MC-27803	930
SC5011	N8120511	Elan Sampler Cone O-rings, pkt 5	15
SC5012	N8125012	Elan Skimmer Cone O-rings, pkt 5	15
SC6000	WE01-7142	Insertion / Removal Tool	210
2085-2206-05		Diamond paste for polishing and cleaning cones, 3-micron, 5-ml syringe	55
2085-2206-18		Diamond paste for polishing and cleaning cones, 3-micron, 18-ml syringe	93
EPOND N°	OEM P/No.	Description, tubes	CHF
See web site		All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12	40
See web site		Flared' or 'Cone-terminated' peri-pump tubes , 2 tag Tygon P-P tubes, pkt 10	40
		Pharmed 2-stop peri-pump tubes, packet of 12	110
MPP-038-PVC	N0777042	Tygon flared 2-stop tubing for MP2 peri-pump, ID=0,38mm, org-green	44
MPP-130-PHR	N0777444	Santoprene 2-stop tubing for MP2 peri-pump, ID=1,30mm, grey-grey	82
PVCSD-30	N0680375	PVC 2-stop pump tubing, 0.030 (0.762mm) ID, Black/Black, 12/pk.	28
PVCSD-30-F	00473550	PVC 2-stop pump tubing, 0.030 (0.762mm) ID, Black/Black, Flared,12/pk.	48
SolFSD-30	N0680376	SolventFlex 2-stop pump tubing, 0.030 (0.762mm) ID, Black/Black, 12/Pk.	36
SolFSD-30-F	N8122012	SolventFlex 2-stop pump tubing, 0.030 (0.762mm) ID, Black/Black, flared, 12/Pk.	36
ERT-040		EPOND 'RubeTubes', Tygon-core sample tube connectors, packet of 40	165
Autosamplers and ICP / AAS Standard and Reference Solutions			
SC-8000-XXXX		Complete range of ESI SC-series autosamplers - please ask for details	POA
ICP-standards		Certified ICP / AAS Standard solutions - http://www.e-pond.biz/shop	
		<i>reference solutions made to measure - on request</i>	
www. e-pond. biz			