

## Standard 2-Stop Tubing for Peristaltic Pumps (Packs of 12)



i.d.	Stop Colors		Non-Flared		Flared*
			PVC 2-stop	Santoprene 2-stop	PVC 2-stop
0.13mm	Orange	Black	ES-4398-1013		ES-4397-1013
0.19mm	Orange	Red	ES-4398-1019		ES-4397-1019
0.27mm	Orange	Blue	ES-4398-1027	ES-4398-2027	ES-4397-1027
0.38mm	Orange	Green	ES-4398-1038	ES-4398-2038	ES-4397-1038
0.44mm	Green	Yellow	ES-4398-1044		ES-4397-1044
0.51mm	Orange	Yellow	ES-4398-1051	ES-4398--2051	ES-4397-1051
0.57mm	White	Yellow	ES-4398-1057		ES-4397-1057
0.64mm	Orange	White	ES-4398-1064	ES-4398-2064	ES-4397-1064
0.76mm	Black	Black	ES-4398-1076	ES-4398-2076	ES-4397-1076
0.89mm	Orange	Orange	ES-4398-1089	ES-4398-2089	ES-4397-1089
0.95mm	White	Black	ES-4398-1095		ES-4397-1095
1.02mm	White	White	ES-4398-1102	ES-4398-2102	ES-4397-1102
1.09mm	White	Red	ES-4398-1109		ES-4397-1109
1.14mm	Red	Red	ES-4398-1114	ES-4398-2114	ES-4397-1114
1.22mm	Red	Grey	ES-4398-1122		ES-4397-1122
1.30mm	Grey	Grey	ES-4398-1130	ES-4398-2130	ES-4397-1130
1.42mm	Yellow	Yellow	ES-4398-1142	ES-4398-2142	ES-4397-1142
1.52mm	Yellow	Blue	ES-4398-1152	ES-4398-2152	ES-4397-1152
1.65mm	Blue	Blue	ES-4398-1165	ES-4398-2165	
1.75mm	Blue	Green	ES-4398-1175		
1.85mm	Green	Green	ES-4398-1185	ES-4398-2185	
2.06mm	Purple	Purple	ES-4398-1206	ES-4398-2206	
2.20mm	Purple	Black	ES-4398-1220	ES-4398-2220	
2.54mm	Purple	Orange	ES-4398-1254	ES-4398-2254	
2.79mm	Purple	White	ES-4398-1279	ES-4398-2279	
3.17mm	Black	White	ES-4398-1317		

\*For easy insertion of Teflon® capillaries **MPP-idd** idd = internal diameter

